
DWX-42W

Setup Guide

Read this first

This guide describes the setup tasks and important conditions about the installation site that must be met in order to enable use of this machine.

Thank you very much for purchasing this product.

- To ensure correct and safe usage with a full understanding of this product's performance, please be sure to read through this manual completely and store it in a safe location.
- Unauthorized copying or transferal, in whole or in part, of this manual is prohibited.
- The contents of this operation manual and the specifications of this product are subject to change without notice.
- The operation manual and the product have been prepared and tested as much as possible. If you find any misprint or error, please inform Roland DG Corporation.
- DGSHAPE Corporation assumes no responsibility for any direct or indirect loss or damage which may occur through use of this product, regardless of any failure to perform on the part of this product.
- DGSHAPE Corporation assumes no responsibility for any direct or indirect loss or damage which may occur with respect to any article made using this product.

Applicable EMC Standards: EN 61326-1 (Class A), EN 55011 (Class A, Group 1)

This product is a Class A device.

It may generate radio wave interference when used in a residential environment. In this situation, the user may be required to implement appropriate countermeasures.

This product is a Group 1 device.

Radio frequency energies in the form of electromagnetic radiation, induction, and/or electrostatic bonding are not used nor deliberately emitted when processing material or during examination/analysis.

Company names and product names are trademarks or registered trademarks of their respective holders.

<http://www.dgshape.com/>

Copyright © 2018 DGSHAPE Corporation

For the USA

**FEDERAL COMMUNICATIONS COMMISSION RADIO
FREQUENCY INTERFERENCE STATEMENT**

Responsible Party: Roland DGA Corporation
Address: 15363 Barranca Parkway Irvine, CA 92618 U.S.A.
Telephone: 949-727-2100
Type of Equipment: Engraving Machine
Model Name: DWX-42W

NOTE:

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC CAUTION

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Use only I/O cables that have been designed and manufactured specifically for this device.

For California

WARNING: This product can expose you to chemicals including lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

For Canada

CAN ICES-3 (A)/NMB-3(A)

For EU Countries

Manufacturer:
ROLAND DG CORPORATION
1-6-4 Shinmiyakoda, Kita-ku, Hamamatsu-
shi, Shizuoka-ken, 431-2103 JAPAN

The importer in the EU:
Roland DG Europe Holdings B.V.
Prof. J.H. Bavincklaan 2, 1183 AT, Amstelveen, The
Netherlands

NOTICE

Grounding Instructions

In the event of a malfunction or breakdown, grounding provides a path of least resistance for electric current to reduce the risk of electric shock. This tool is equipped with an electric cord having an equipment-grounding conductor and a grounding plug. The plug must be plugged into a matching outlet that is properly installed and grounded in accordance with all local codes and ordinances.

Do not modify the plug provided - if it will not fit the outlet, have the proper outlet installed by a qualified electrician.

Improper connection of the equipment-grounding conductor can result in a risk of electric shock. The conductor with insulation having an outer surface that is green with or without yellow stripes is the equipment-grounding conductor. If repair or replacement of the electric cord or plug is necessary, do not connect the equipment-grounding conductor to a live terminal.

Check with a qualified electrician or service personnel if the grounding instructions are not completely understood, or if in doubt as to whether the tool is properly grounded.

Use only 3-wire extension cords that have 3-prong grounding plugs and 3-pole receptacles that accept the tool's plug.

Repair or replace damaged or worn cord immediately.

Operating Instructions

KEEP GUARDS IN PLACE and in working order.

MAKE WORKSHOP KID PROOF with padlocks, master switches, or by removing starter keys.

DON'T OVERREACH. Keep proper footing and balance at all times.

DISCONNECT TOOLS before servicing; when changing accessories, such as blades, bits, cutters, and the like.

REDUCE THE RISK OF UNINTENTIONAL STARTING. Make sure switch is in off position before plugging in.

NEVER STAND ON TOOL. Serious injury could occur if the tool is tipped or if the cutting tool is unintentionally contacted.

CHECK DAMAGED PARTS. Before further use of the tool, a guard or other part that is damaged should be carefully checked to determine that it will operate properly and perform its intended function - check for alignment of moving parts, binding of moving parts, breakage of parts, mounting, and any other conditions that may affect its operation. A guard or other part that is damaged should be properly repaired or replaced.

NEVER LEAVE TOOL RUNNING UNATTENDED. TURN POWER OFF. Don't leave tool until it comes to a complete stop.

USE POWER SUPPLY CORD WHICH IS ATTACHED WITH PRODUCTS. do not use other power supply cord.

Contents

Contents	2
To Ensure Safe Use	4
Pour utiliser en toute sécurité	12
Important Notes on Handling and Use	20
1. Getting Started	21
About This Machine	22
Features of This Machine	22
Documentation Included with the Machine	22
Checking the Part Names and Included Items	23
Part Names and Functions	23
Checking the Included Items	27
Detection Pin for Correction	28
Items to Prepare Yourself	29
Items to Prepare Yourself	29
Water (Soft or Purified Water)	29
Source of Compressed Air (Compressor)	29
Air Hose	30
Fluid Receptacle	30
Additive (ZAW-1000D)	31
Chelating Reagent (ZCH-100D)	31
2. Installation and Setup	33
1. Installation	34
Step 1 : Deciding on an Installation Site	34
Step 2 : Removing the Retainers	36
Step 3 : Attaching the Regulator	37
Step 4 : Attaching the Garbage Receptacle	42
2. Connecting to the Computer	43
Procedure Overview	43
Step 1 : Installing the Software	44
Step 2 : Connecting the USB Cable	47
Step 3 : Switching the Power On	48
Step 4 : Selecting the Machine to Connect To	50
Step 5 : Determining the LAN Address	51
Step 6 : Setting the Computer's Network	52
Step 7 : Setting the Machine's Network	54
Step 8 : Disconnecting the USB Cable	57
Step 9 : Connecting the LAN Cable	57
Step 10 : Selecting the Machine to Connect To	59
3. Performing the Initial Work	61
Initial Work	61
Step 1 : Filling the Coolant Tank	61
Step 2 : Automatic Correction	68
Step 3 : Spindle Run-In (Long)	73

3. Appendix.....	77
Specifications.....	78
Power Rating and Serial Number Location	78
Coordinates	78
Dimensional Drawings	79
Specifications.....	81

To Ensure Safe Use

Improper handling or operation of this machine may result in injury or damage to property. Points which must be observed to prevent such injury or damage are described as follows.

About WARNING and CAUTION Notices

 WARNING	<p>Used for instructions intended to alert the user to the risk of death or severe injury should the unit be used improperly.</p>
 CAUTION	<p>Used for instructions intended to alert the user to the risk of injury or material damage should the unit be used improperly.</p> <p>* Material damage refers to damage or other adverse effects caused with respect to the home and all its furnishings, as well as to domestic animals or pets.</p>

About the Symbols

	<p>The  symbol alerts the user to important instructions or warnings. The specific meaning of the symbol is determined by the design contained within the  symbol. The symbol at left means “danger of electrical shock.”</p>
	<p>The  symbol alerts the user to items that must never be carried out (are forbidden). The specific meaning of the symbol is determined by the design contained within the  symbol. The symbol at left means the unit must never be disassembled.</p>
	<p>The  symbol alerts the user to things that must be carried out. The specific meaning of the symbol is determined by the design contained within the  symbol. The symbol at left means the power-cord plug must be unplugged from the outlet.</p>

 **This is a heavy machine**

CAUTION

Unloading and emplacement are operations that must be performed by 2 persons or more.

Tasks that require undue effort when performed by a small number of persons may result in physical injury. Also, if dropped, such items may cause injury.

CAUTION

Install in a level and stable location.
Failure to do so may result in falling of the machine, leading to injury.

 **Incorrect operation may cause injury**

 **WARNING**

Always unplug the power cord when attaching or removing parts and optional parts and when performing maintenance that does not require the machine to be connected to a power source.

Attempting such operations while the machine is connected to a power source may result in injury or electrical shock.

Keep children away from the machine.

The machine includes areas and components that pose a hazard to children and may result in injury, blindness, choking, or other serious accident.

Never attempt to disassemble, repair, or modify the machine.

Doing so may result in fire, electrical shock, or injury. Entrust repairs to a trained service technician.

Be sure to follow the operation procedures described in the user's manual. Never allow anyone unfamiliar with the usage or handling of the machine to touch it.

Incorrect usage or handling may lead to an accident.

Never operate the machine while tired or after ingesting alcohol or any medication.

Operation requires unimpaired judgment. Operating with impaired judgment may result in an accident.

Conduct operations in a clean, brightly lit location.

Working in a location that is dark or cluttered may lead to an accident, such as becoming caught in the machine as the result of an inadvertent stumble.

Never climb or lean on the machine.

The machine is not made to support a person. Climbing or leaning on the machine may dislodge components and cause a slip or fall, resulting in injury.

 **WARNING**

Be careful around the milling tool.

The milling tool is sharp. Broken milling tools are also dangerous. To avoid injury, exercise caution.

Never use a milling tool that has become dull. Perform frequent maintenance to keep the machine in good working order.

Inappropriate usage may result in fire or injury.

Never operate the machine if the cover is cracked or broken.

Doing so may result in injury. If the cover is cracked or broken, contact your authorized Roland DG Corporation dealer.

Be careful of the pointed portion inside the front cover.

There is a pointed portion inside the front cover. Exercise caution when cleaning.

 **CAUTION**

Use a dry cloth to clean the inside of the equipment.

Failure to do so may cause the components inside the equipment to degrade, which can lead to injury.

Wipe up any water around the main unit.

Failure to do so may cause a slip or fall, resulting in injury.

 **Milling waste and workpieces are flammable and toxic.**

 **WARNING**

Never attempt to mill magnesium or any other such flammable material.
Fire may occur during milling.

Keep open flame away from the work area.
Milling waste may ignite. Powdered material is extremely flammable, and even metal material may catch fire.

Never use a vacuum cleaner to clean up milling waste.
Picking up fine cuttings using an ordinary vacuum cleaner may result in fire or explosion.

 **CAUTION**

Wear dust goggles and a mask. Wash away any milling waste remaining on the hands.

Accidentally swallowing or inhaling milling waste may be hazardous to your health.

 **Danger of pinching, entanglement, and burns**

 **WARNING**

Securely fasten the milling tool and workpiece in place. After securing in place, make sure no spanners or other articles have been left behind inadvertently.
Otherwise, such articles may be thrown from the machine with force, posing a risk of injury.

Do not touch the spindle unit or the surrounding areas immediately after milling has ended.
Doing so may result in burns.

 **CAUTION**

Never attempt operation while wearing a necktie, necklace, or loose clothing. Bind long hair securely. When performing operations while wearing gloves, exercise sufficient caution to prevent entanglement in the machine.
Such items may become caught in the machine, resulting in injury.


 **CAUTION**


Do not use hot water for the coolant.
Using hot water may cause burns.

When opening or closing the front cover, hold the part indicated in the following figure with both hands.
Failing to do so may result in your fingers being pinched, leading to injury.

When raising or lowering the coolant tank, hold the parts indicated in the following figure.
Failing to do so may result in your fingers being pinched, leading to injury.

Danger of electrical short, electrical shock, or fire

WARNING

Connect to an electrical outlet that complies with this machine's ratings (for voltage, frequency, and current). Incorrect voltage or insufficient current may cause fire or electrical shock.

If sparking, smoke, burning odor, unusual sound, or abnormal operation occur, immediately unplug the power cord. Never use if any component is damaged. Continuing to use the machine may result in fire, electrical shock, or injury. Contact your authorized Roland DG Corporation dealer.

Never use outside or in any location where exposure to water or high humidity may occur. Do not touch the power cord, plug, or electrical outlet with wet hands. Doing so may result in fire or electrical shock.

Never allow any foreign object to get inside. Never expose the machine to liquid spills. Inserting objects such as coins or matches or allowing beverages to be spilled into the ventilation ports may result in fire or electrical shock. If anything gets inside, immediately disconnect the power cord and contact your authorized Roland DG Corporation dealer.

Never place any flammable object nearby. Never use a combustible aerosol spray nearby. Never use in a location where gases can accumulate. Combustion or explosion could occur.

Handle the power cord, plug, and electrical outlet correctly and with care. Never use any article that is damaged. Using a damaged article may result in fire or electrical shock.

When using an extension cord or power strip, use one that adequately satisfies the machine's ratings (for voltage, frequency, and current). Use of multiple electrical loads on a single electrical outlet or of a lengthy extension cord may cause fire.

WARNING

Position the machine so that the power plug is within immediate reach at all times.

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For accessories (optional and consumable items, power cord, and the like), use only genuine articles compatible with this machine.

Incompatible items may lead to an accident.

Never use the machine for any purpose for which it is not intended, or use the machine in an excessive manner that exceeds its capacity. Doing so may result in injury or fire.

Do not use the supplied power cord for other products.

If the machine will not be used for a long period, disconnect the power cord. This prevents an accident caused by unexpected electric leakage or unintentional start of the machine.

Connect to ground. This can prevent fire or electrical shock due to current leakage in the event of a malfunction.

Never use cutting oil. This machine is not designed for the flow of cutting oil. Oil may get inside the machine and cause fire or electrical shock.

Never use a pneumatic blower. This machine is not compatible with a pneumatic blower. Milling waste may get inside the machine and cause fire or electrical shock.

⚠ WARNING

Handle the coolant filter carefully.
Dropping the filter or grasping it too strongly may damage the filter. If the filtration function does not work due to damage, the line may be clogged with milling waste, which may lead to coolant leaks and fire and electrical shock due to current leakage.


Use a dry cloth when cleaning silicone resin parts, and be careful to prevent the silicone resin from being damaged.
The silicone resin being damaged may lead to electric leakage.

⚠ WARNING


Do not place any electronic devices in the vicinity of this machine.
Because coolant flows within this machine, water may be sprayed on objects in the vicinity of the machine when its front cover is opened. To prevent malfunctions, do not place any electronic devices in the vicinity of this machine.


⚠ Important notes about the power cord, plug, and electrical outlet


Never place any object on top or subject to damage.

Never bend or twist with undue force.


Never pull with undue force.


Never bundle, bind, or roll up.

Never allow to get wet.

Never make hot.

Dust may cause fire.

Danger of components flying off or becoming damaged

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Never use in a location exposed to cutting oil, solvents, chemicals, or other such substances.

Components may deteriorate or rupture under the pressure of compressed air.

Never use in a location exposed to direct sunlight.

Components may deteriorate or rupture under the pressure of compressed air.

Keep the temperature of the installation area within the specified range. Never place a stove or heater nearby.

Components may deteriorate or rupture under the pressure of compressed air.

Connect the air hose securely, so that it will not come loose.

A pressurized hose that comes loose may whip about uncontrolled, posing a hazard. Connect securely.

Before attempting to attach or detach the air hose or optional items or before attempting to perform cleaning or maintenance that does not require the machine to be connected to a power source, stop the supply of compressed air and allow the pressure to escape.

Allowing the equipment to remain under pressure poses a hazard of flying components.

Use an air hose of the specified diameter and having adequate pressure resistance.

Otherwise the article may come loose or rupture.

Keep the supply of compressed air within the specified pressure range.

Exceeding the specified pressure may result in a rupture or other major accident.

WARNING

Never strike or subject the object to impact.

Components may suffer damage or rupture under the pressure of compressed air.

When the machine will be out of use for a prolonged period, stop supply of compressed air and bleed off the air pressure.

This can prevent an accident.

Never damage the air hose, or bend or twist it with undue force. Never use an item that has deteriorated.

A damaged air hose may rupture.


Be sure to bleed off the air pressure before removing the regulator bowl.

Failure to do so may result in a rupture or components flying off.



Clean the regulator bowl using a neutral detergent. Never use solvents such as gasoline, alcohol, or thinner.

Using solvents may degrade the bowl and may result in a rupture.


⚠ Important notes about the end (cut edge) of the air hose


Failure to observe these cautions may result in air leakage or easy dislodgement of the air hose. Also, when a hose has been detached, cut off the end before reattaching the hose.

Ensure that it is not crushed.

Make the cut edge straight.

Keep free of damage and soiling.

⚠ Additives and coolant are toxic

⚠ WARNING

Never drink or sniff the additives or coolant or allow them to come in contact with the eyes or skin.
Doing so may be hazardous to your health.

⚠ WARNING

After milling, wash away any products with purified water or the like.
After milling, there will be coolant on the product. Using the product as is may cause an inflammation or the like.

In the event of ingestion or physical distress

- In the event of contact with the eyes, immediately flush with water for at least 15 minutes. If eye irritation continues, seek treatment by a physician.
 - In the event of contact with skin, immediately wash with soap. If irritation or inflammation occurs, seek treatment by a physician.
 - In the event of ingestion, do not induce vomiting, and immediately seek treatment by a physician. Forcibly inducing vomiting may lead to danger of choking.
 - If odor leads to physical distress, move to a well-ventilated location and rest quietly. If dizziness or nausea persists, seek treatment by a physician.
-
- Refer to the appropriate safety data sheet (SDS) for the chemical substances used in the grease and the safety related to those substances.
 - Refer to the appropriate safety data sheet (SDS) for the chemical substances used in the additives and the safety related to those substances.

 **Handle the machine with care**

 **CAUTION**

Do not overfill or tilt the coolant tank.
The fluid inlet on the coolant tank is open.
Overfilling or tilting will cause fluid to spill out.


 **CAUTION**


Be sure to put the specified additive in the coolant tank.
Additives are effective in reducing coolant deterioration and raising the milling efficiency in order to maintain product performance. If additives are not used, the coolant may generate an unpleasant odor.

 **Warning labels**

Warning labels are affixed to make areas of danger immediately clear. The meanings of these labels are as follows. Be sure to heed their warnings. Also, never remove the labels or allow them to become obscured.

Caution: High Voltage
Removing the cover may result in high-voltage electrical shock.

Caution: Sharp Tool
The milling tool is sharp. Inadvertent contact may cause injury.

Never use a pneumatic blower.
This machine is not compatible with a pneumatic blower. Milling waste may get inside the machine and cause fire or electrical shock.

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Important Notes on Handling and Use

This machine is a precision device. To ensure the full performance of this machine, be sure to observe the following important points. Failure to observe these points may not only result in loss of performance but may also cause malfunction or breakdown.

This machine is a precision device.

- Handle this machine carefully and never subject it to impacts or excessive force.
- Carefully clean away milling waste.
- Use this machine within the range of its specifications.
- Never attempt to move the spindle unit or rotary axis by hand with undue force.
- Never needlessly touch anywhere inside the machine except for locations specified in this manual.

Install this machine in a suitable location.

- Install this machine in a location that meets the specified conditions for temperature, humidity, and the like.
- Install this machine in a quiet, stable location offering good operating conditions.
- Never install this machine outside.
- Never install this machine in an environment where silicone substances (oil, grease, spray, etc.) are present. Doing so may cause poor switch contact.

Leave the machine and the compressor on at all times.

- Turning the machine off leads to moisture collecting inside of the machine, which may lead to malfunction.
- If the compressor is turned off, the supply of air to the machine's spindle unit will be stopped, which may lead to malfunction.

This machine becomes hot.

- Never cover the ventilation holes with cloth, tape, or anything else.
- Install this machine in a well-ventilated location.

This machine is a dedicated cutting machine for glass ceramics and composite resins.

- Never cut any material other than glass ceramics and composite resins.

1. Getting Started

- About This Machine 22
 - Features of This Machine 22
 - Documentation Included with the Machine..... 22
- Checking the Part Names and Included Items..... 23
 - Part Names and Functions 23
 - Checking the Included Items 27
 - Detection Pin for Correction 28
- Items to Prepare Yourself..... 29
 - Items to Prepare Yourself 29
 - Water (Soft or Purified Water) 29
 - Source of Compressed Air (Compressor)..... 29
 - Air Hose 30
 - Fluid Receptacle..... 30
 - Additive (ZAW-1000D)..... 31
 - Chelating Reagent (ZCH-100D) 31

About This Machine

Features of This Machine

This machine is a dental cutting and milling machine that uses glass ceramics and composite resins to make dental prostheses, including crowns, bridges, inlays, onlays, and veneers.

Installing and setting up the exclusive software on your computer and connecting your computer to this machine enables you to create high-quality dental prostheses.

For the latest information regarding this machine (including manuals), see the DGSHAPE Corporation website (<http://www.dgshape.com>).

Documentation Included with the Machine

The following documentation is included with the machine.

- **DWX-42W Setup Guide (this document)**

This manual describes how to install the machine and provides information on such matters as how to install the included software. Be sure to read this first.

* You can find the PDF of the manual on the included DGSHAPE Software Package CD. The contents of the printed document and the PDF are the same.

- **DWX-42W User's Manual (electronic-format manual)**

This manual explains how to operate the software and how to perform milling using this machine. The manual is in electronic format, and no printed document is included. You can find the manual on the included DGSHAPE Software Package CD.

☞ P. 61 "Step 1 : Filling the Coolant Tank"


- **Daily Maintenance Card ("Before/After Daily Operations")**

This document provides an overview of the operations that must be performed before and after daily operations. Be sure to perform the operations on this card before starting and after finishing machine operations each day.


Checking the Part Names and Included Items

Part Names and Functions

Front

Back


MEMO


The regulator can also be installed on the side of the machine.

☞ P. 37 "Step 3 : Attaching the Regulator"

Inside the Front Cover


Rotary axis unit		
④	Dummy pin stocker	This is a stocker dedicated for the dummy pin. Do not set a milling bur in this stocker. Also, do not remove the dummy pin that is returned to this stocker.
⑤	Milling bur sensor	This sensor is used to adjust the milling bur height. If the surface is dirty, the sensor's accuracy will decrease, so keep the surface clean.
⑥	Stocker	The milling burs are set here.
⑦	Mounting part	The workpiece is attached here. You can attach up to three workpieces.
⑧	ATC magazine	Up to six milling burs can be set in this ATC-compliant magazine. ATC stands for Auto Tool Changer. This is a function for automatically replacing the milling bur (tool) to match the milling data.

Spindle head			
	①	Coolant nozzle	The coolant is sprayed from this nozzle. Clean it once a week. ☞ "User's Manual"
	②	Collet	This tool is used to secure the milling bur. It needs to be maintained after operations are finished each day. ☞ "User's Manual"
	③	Dummy pin	This is a milling bur dummy. It is attached to the spindle upon shipment. A malfunction may occur if the collet does not grasp anything, so have it grasp the dummy pin when it is not grasping a milling bur.


* In the manuals for this machine (including this document), the mechanisms around the spindle unit are called the "spindle head." Also, the rotary-axis area inside the spindle unit is called the "spindle."

Inside the Bottom Cover

















①	Coolant tank	This is the tank that holds the coolant. The coolant must be replaced and the filters must be cleaned periodically.
②	Storage	Milling burs and other items can be stored here. It has two levels. The upper level is exactly large enough to fit the included maintenance kit case.

Side

Checking the Included Items

The following items are included with the machine. Make sure they are all present and accounted for.

 <p>Power cord</p>	 <p>USB cable (1)</p>	 <p>Drain hoses (2)</p>	 <p>Dedicated milling burs (4)</p>	
 <p>Spare tool retaining screws (3)</p>	 <p>Rotary axis correction jig (1)</p>	 <p>ATC magazine correction jig (1)/ A screw (1)</p>	 <p>Detection pin for correction (1)</p>	
 <p>Measuring cup (1)</p>	 <p>Coolant filters (2)</p>	 <p>Regulator (with stay)/ B screws (2)</p>	 <p>Waste receptacle</p>	
 <p>Air hose (1)</p>	 <p>Milling bur removal jig (1)</p>	 <p>Hexagonal screwdriver (1)</p>	 <p>Cloth for care and maintenance</p>	 <p>Setup Guide (this document; 1)/ Daily Maintenance Card (1)</p>
 <p>Maintenance kit</p> <ul style="list-style-type: none"> ① Case (1) ② Grease (1) ③ Small case (1) ④ Collet brush (1) ⑤ Nozzle cleaners (2) ⑥ Taper cleaner (3) ⑦ Collet tap (1) ⑧ Collet replacement jig (1) 			 <p>DGSHAPE Software Package CD (1)</p>	

Detection Pin for Correction

- The detection pin for correction is a maintenance tool used to correct the position of the rotary axis.
- Handle the “detection pin for correction” and the “dummy pin” separately to avoid confusing them.

The “dummy pin” is used as a “dummy” milling bur. A malfunction may occur if the collet does not grasp anything, so have it grasp the dummy pin when it is not grasping a milling bur. The dummy pin is attached to the spindle upon shipment.

The dummy pin and the detection pin for correction have the same shape, but do not confuse them. If the detection pin for correction is used as the dummy pin even once, the detection pin for correction can no longer be used to provide proper correction. If for some reason the detection pin for correction is used as the dummy pin, a new detection pin for correction will be necessary. Contact your authorized DGSHAPE Corporation dealer or access our website (<http://www.dgshape.com/>).

Items to Prepare Yourself

Items to Prepare Yourself

- Water (soft or purified water)
- Source of compressed air (compressor)
- Air hose
- Fluid receptacle
- Additive (ZAW-1000D)
- Chelating reagent (ZCH-100D)

Explanations of each item are given below.

Water (Soft or Purified Water)

This is required in order to make the coolant. Use soft or purified water. Using hard water may have a negative effect on the service life of the milling bur and on the quality of the product.

Source of Compressed Air (Compressor)

This machine requires compressed air. You will need to prepare a compressor separately.

WARNING

The pressure of the compressed air must be 1.0 MPa or less.

Anything higher may result in a serious accident such as a rupture.

WARNING

Ensure that the supplied compressed air is not contaminated by water, oil, chemicals, or foreign objects.

The components may deteriorate or rupture, or the contaminants may be scattered, posing a hazard.

WARNING

Do not supply the compressed air until the air hose is connected securely.

Failure to do so may cause an accident.

The compressor must meet the following conditions.

- Control pressure: 0.3 MPa or more
- Air capacity: 39 L/min or more
- Oil-free type
- Dryer-equipped


Air Hose

You will need to separately prepare an air hose to connect the compressor to the regulator.

*** Use the included air hose to connect the regulator to the machine.**

The air hose must meet the following conditions.

- Hose outer diameter: 6 mm (polyurethane resin tube)


Note

Be sure to use the included regulator when supplying compressed air.

Fluid Receptacle

Because fluid accumulated in the regulator will drain out little by little, prepare a fluid receptacle such as a tray.

☞ P. 37 "Step 3 : Attaching the Regulator"

Additive (ZAW-1000D)

To make the coolant, you will need to prepare the designated additive separately. To purchase additives, contact your authorized DGSHAPE Corporation dealer or access our website (<http://www.dgshape.com/>).

Note

- Store additives in a cool, dark place.
- Due to the characteristics of the internal components, the color may change and deposits of separated components may form, but these do not indicate any problems with the use of the additives.
- If component deposits form, lightly shake the container to mix the components before use.
 - * **The expiration date is one year from the date of manufacture (written on the unit packaging box).**

Chelating Reagent (ZCH-100D)

The ZCH-100D (hereinafter referred to as the “chelating reagent”) works to make calcium and other substances contained in water inactive. Including it in this product makes it easier to clean off milling powder that adheres to the machine and the cutting tool after milling. To purchase chelating reagent, contact your authorized DGSHAPE Corporation dealer or access our website (<http://www.dgshape.com/>).

* **The expiration date is one year from the date of manufacture (written on the unit packaging box).**

In the event of additive or chelating reagent ingestion or physical distress

- In the event of contact with the eyes, immediately flush with water for at least 15 minutes. If eye irritation continues, seek treatment by a physician.
- In the event of contact with skin, immediately wash with soap. If irritation or inflammation occurs, seek treatment by a physician.
- In the event of ingestion, do not induce vomiting, and immediately seek treatment by a physician. Forcibly inducing vomiting may lead to danger of choking.
- If odor leads to physical distress, move to a well-ventilated location and rest quietly. If dizziness or nausea persists, seek treatment by a physician.
 - * **Refer to the appropriate safety data sheet (SDS) for the chemical substances used in the additive and chelating reagent and the safety related to those substances.**

2. Installation and Setup

1. Installation.....	34
Step 1 : Deciding on an Installation Site.....	34
Step 2 : Removing the Retainers.....	36
Step 3 : Attaching the Regulator.....	37
Step 4 : Attaching the Garbage Receptacle.....	42
2. Connecting to the Computer.....	43
Procedure Overview.....	43
Step 1 : Installing the Software.....	44
Step 2 : Connecting the USB Cable.....	47
Step 3 : Switching the Power On.....	48
Step 4 : Selecting the Machine to Connect To.....	50
Step 5 : Determining the LAN Address.....	51
Step 6 : Setting the Computer's Network.....	52
Step 7 : Setting the Machine's Network.....	54
Step 8 : Disconnecting the USB Cable.....	57
Step 9 : Connecting the LAN Cable.....	57
Step 10 : Selecting the Machine to Connect To.....	59
3. Performing the Initial Work.....	61
Initial Work.....	61
Step 1 : Filling the Coolant Tank.....	61
Step 2 : Automatic Correction.....	68
Step 3 : Spindle Run-In (Long).....	73

1. Installation

Step 1 : Deciding on an Installation Site

Placement and Installation


The weight of the machine alone is 61 kg (135 lb.). Unload and place the machine with care.

⚠ WARNING

Unloading and placement are operations that must be performed by 2 persons or more. Tasks that require undue effort when performed by a small number of people may result in physical injury. Also, if dropped, such items may cause injury.

⚠ CAUTION

When lifting the machine, hold the locations shown in the figure.

Installation Environment

⚠ WARNING

Install in a level and stable location.

Failure to do so may result in falling of the machine, leading to injury.

⚠ WARNING

Never install the machine in a location exposed to open flame.

Milling waste may ignite. Powdered material is extremely flammable, and even metal material may catch fire.

⚠ WARNING

Never install the machine close to any flammable object or in a gas-filled location.

Combustion or explosion may occur.

⚠ WARNING

Never install this machine outside or in any location where exposure to water or high humidity may occur.

Current leakage may cause electrical shock or fire.

⚠ WARNING


Position so that the power plug is within immediate reach at all times.

This is to enable quick disconnection of the power plug in the event of an emergency. Install the machine next to an electrical outlet. Also, provide enough empty space to allow immediate access to the electrical outlet.

- Never install the machine in a location subject to wide fluctuations in temperature or humidity.
- Never install the machine in a location subject to shaking or vibration.
- Never install the machine in a dusty or dirty location.
- Never install the machine in a location exposed to direct sunlight or near air-conditioning or heating equipment.
- Never install the machine in a location exposed to considerable electrical or magnetic noise or other forms of electromagnetic energy.
- Never install this machine in an environment where silicone substances (oil, grease, spray, etc.) are present. Doing so may cause poor switch contact.


Installation space

Ensure that at least the following amount of space is available.

Ensure that enough space is available to allow the front cover to be opened.

Ensure that space is available in front of and behind the machine to enable the opening of the bottom cover.

Installation Location Height

- **The installation location should be above the work floor by 0.6 m (23.6 in.) or more.**

This machine is a desktop-type unit. Install in a location where the power button and other controls are easily accessible.

1. Installation


Step 2 : Removing the Retainers

Retainers are attached to protect the machine from vibration during shipment. When installation is complete, remove and save these retainers.

- Be sure to remove the retainers. Any that remain may cause faulty operation or a breakdown when the power is switched on.
- The retainers will be required when moving the machine to a different location. Store them carefully so that they do not get misplaced.
- When removing the retainers, make sure the power cord is not connected.


Procedure

- 1 Open the front cover.


- 2 Remove the retainers.

- 1 Pull the bottom of the retainer (large) toward you.
- 2 Slide the retainer with it facing up to remove it.

- ③ Remove the retainer (small).


- ③ Close the front cover.

Step 3 : Attaching the Regulator

The regulator for this machine is the auto-drain type. Attaching the regulator incorrectly can cause leakage. Be sure to attach the regulator according to the procedure.

⚠ WARNING Before supplying compressed air, make sure that the bowl is securely attached. If the bowl is not properly attached, it may come flying off.

Required Included Items

<p>Regulator (with stay)/ B screws (2)</p>	<p>Hexagonal screwdriver (1)</p>
--	----------------------------------

1. Installation


Procedure

1 Attach the regulator to the machine.


When attaching the regulator to the back of the machine

① Temporarily tighten the B screws in two locations on the back of the machine.

Use the screw holes at the top of the back of the machine.


② Attach the regulator's stay to the B screws from step ①.

③ Secure the regulator's stay.


Fully tighten the B screws attached in step ①.

When attaching the regulator to the side of the machine


① **Remove the regulator from the stay.**

Loosen the screws.

② **Temporarily tighten the B screws in two locations on the side of the machine.**


Use the screw holes at the bottom of the side of the machine. You can attach the regulator on either the left side or right side of the machine.

MEMO

There are screw holes for attaching the regulator at the same height on the left side and right side of the machine.


③ **Attach the regulator to the B screws from step ②.**


1. Installation

④ Secure the regulator.


Fully tighten the B screws attached in step ②.


② Use the air hoses to connect the regulator to the machine and the regulator to the source of compressed air.

To attach: Insert firmly, as far as it will go.

To detach: Press and hold the ring and pull out the hose.

⚠ WARNING

Securely insert the air hose as far as it will go. Lightly tug the hose to make sure it does not come loose.

If it is not inserted securely, it may come loose.


③ Supply compressed air and verify that air does not leak from the connections or any other location.

4 Adjust the air pressure.

⚠ WARNING Turn the air pressure adjustment knob slowly and carefully. Otherwise, the machine may move suddenly, posing a risk of injury.


- ① Pull up on the upper knob (the air pressure adjustment knob).
- ② Slowly turn the upper knob.
Adjust until the meter reads between 0.18 and 0.22 MPa.
- ③ Push down on the upper knob.

Important: Be sure to set the regulator to 0.22 MPa or lower. Anything higher may result in a malfunction.

5 Prepare a fluid receptacle.


Use the drain hoses (included) or a tray.

Step 4 : Attaching the Garbage Receptacle

Procedure

- 1 Open the front cover.
- 2 Attach the garbage receptacle.

2. Connecting to the Computer

Procedure Overview

You can connect the computer and this machine through its USB interface or LAN interface. The procedure varies depending on the connection method. We recommend that you use the overview shown below to understand the overall procedure before starting the connection operations.

USB	LAN	Procedure
<input type="radio"/>	<input type="radio"/>	Step 1: Installing the Software (P. 44)
<input type="radio"/>	<input type="radio"/>	Step 2: Connecting the USB Cable (P. 47)
<input type="radio"/>	<input type="radio"/>	Step 3: Switching the Power On (P. 48)
<input type="radio"/>	<input type="radio"/>	Step 4: Selecting the Machine to Connect To (P. 50)
	<input type="radio"/>	Step 5: Determining the LAN Address (P. 51)
	<input type="radio"/>	Step 6: Setting the Computer's Network (P. 52)
	<input type="radio"/>	Step 7: Setting the Machine's Network (P. 54)
	<input type="radio"/>	Step 8: Disconnecting the USB Cable (P. 57)
	<input type="radio"/>	Step 9: Connecting the LAN Cable (P. 57)
	<input type="radio"/>	Step 10: Selecting the Machine to Connect To (P. 59)

* * When using a LAN connection, we recommend that you use a LAN connection after connecting over a USB according to this procedure. By establishing a USB connection and a LAN connection, you can reduce the risk of operations being stopped due to unexpected network problems.

* * You can connect just over a LAN without using a USB connection. For details on the connection method, see the "User's Manual."

Step 1 : Installing the Software

System Requirements

Operating system (OS)	Windows 7/8.1/10 (32-bit and 64-bit versions)
CPU	Minimum required CPU for the operating system
Memory	Minimum amount of RAM required for the operating system
Optical drive	CD-ROM drive
Video card and monitor	A monitor with at least 256 colors and a resolution of 1024 × 768 or more is recommended.

- This software is a 32-bit application and therefore runs in WOW64 (Windows-On-Windows 64) when running on 64-bit versions of Windows operating systems.
- This software has not been tested when running in the Windows XP mode in Windows 7.

For the latest information, see the DGSHAPE Corporation website (<http://www.dgshape.com/>).

Installable Software

VPanel for DWX	This is the machine control software designed for the DWX series made by DGSHAPE Corporation.* * Some models are not supported. For details on the supported models, contact your authorized Roland DG Corporation dealer.
Dental Driver (Windows driver)	This is a Windows-based driver required for sending data from a computer to the machine.
Setup Guide (electronic-format manual)	This is the PDF of this document. The contents are the same.
User's Manual (electronic-format manual)	This manual explains how to operate the software, how to perform milling using this machine, and how to perform maintenance.

Installation Procedure

You can install the driver, software, and electronic-format manual all at once. You can also install the driver and each piece of software separately.

☞ “User’s Manual” (electronic-format manual) “Installing the Driver Separately”

Important

Be sure to connect the machine to the computer as instructed in the procedure. Failure to follow the correct procedure may make installation impossible.

Procedure

1 Before installation, confirm that the machine and the computer are not connected with the USB cable.

2 Log on to Windows as an “Administrator.”

3 Insert the DGSHAPE Software Package CD into the CD-ROM drive of the computer.


When the automatic playback window appears, click [Run menu.exe]. If a [User Account Control] window appears, click [Allow] or [Yes], and then continue with the installation. The setup menu screen appears automatically.

4 Click [Install].

Install the DWX-42W driver and the various software all at once. Follow the on-screen instructions to proceed with the installation.


5 Click [OK].


2. Connecting to the Computer

- 6 Select [Install], [DGSHAPE DWX-42W] for the model name, and [USB] for the port, and then click [Start].

The driver installation begins. Proceed with the installation by following the displayed information.


- 7 If the following window appears, click [Install].

Installation information for each piece of software will be displayed automatically. Follow the instructions displayed for each installer to continue with the installation.

- 8 Click [OK].

- 9 Click X on the setup menu screen.
- 10 Remove the DGSHAPE Software Package CD from the computer.

How to View the Electronic-Format Manual

Windows 10 and 7


From the [Start] menu, click [All Apps] (or [All Programs]), then [DGSHAPE DWX-42W]. Then click [User's Manual].

Windows 8.1

On the [Start] screen, click . On the [Apps] screen, click the [User's Manual] icon under [DGSHAPE DWX-42W].

Step 2 : Connecting the USB Cable

- Use the included USB cable.
- Do not use a USB hub. Connection may not be possible.
- When using a LAN connection, we recommend that you use a LAN connection after connecting over a USB. By establishing a USB connection and a LAN connection, you can reduce the risk of operations being stopped due to unexpected network problems.
- Fix the USB cable in place with the cable clamp.

Important: When connecting multiple machines

When using USB cables to connect two or more machines to one computer, be sure to refer to the User's Manual. To connect the second and subsequent machines, you have to change the machine's ID. Connecting before changing the ID may render the machine unusable.

Step 3 : Switching the Power On

Procedure

1 Connect the power cord.

⚠ WARNING Connect to an electrical outlet that complies with this machine's ratings (for voltage, frequency, and current).

Incorrect voltage or insufficient current may cause fire or electrical shock.

⚠ WARNING Handle the power cord, plug, and electrical outlet correctly and with care. Never use any article that is damaged.


Using a damaged article may result in fire or electrical shock.

⚠ WARNING When using an extension cord or power strip, use one that adequately satisfies the machine's ratings (for voltage, frequency, and current).

Use of multiple electrical loads on a single electrical outlet or use of a lengthy extension cord may cause fire.


⚠ WARNING Connect to ground. This can prevent fire or electrical shock due to current leakage in the event of a malfunction.

⚠ WARNING Connect this machine to a power outlet. Never connect this machine directly to a power distribution panel or other such fixed wiring equipment. Doing so increases the risk of fire or electrical shock.

2 Turn on the power.

The machine starts the initial operations. The Windows driver will be installed automatically on the computer connected with the USB cable.

Important

- Do not bind the power cord with the cable clamp.
- Wire the power cord so that it does not come into contact with the USB cable. The power cord coming into contact with any other cables will result in communication malfunctions due to electrical noise.

Step 4 : Selecting the Machine to Connect To

1 Start VPanel.

Windows 10 and 7


From the [Start] menu, click [All apps] (or [All Programs]), [VPanel for DWX], and then click [VPanel for DWX].

Windows 8.1

Click  on the [Start] screen. From the [Apps] screen, click the [VPanel for DWX] icon under [VPanel for DWX].

2 From the [Usable machine list], click the image of the machine that you want to operate.

The clicked machine is displayed under [Machine to operate] and is displayed dimly in the [Usable machine list]. If you click the same machine again in the [Usable machine list], the machine disappears from the [Machine to operate] field.

3 Click [OK].

- * For a USB connection, this completes the connection to the computer.
- * For a LAN connection, proceed to P. 51 "Step 5 : Determining the LAN Address".

Step 5 : Determining the LAN Address

Important: Windows firewall function

During LAN communication, this machine uses TCP ports 9100, 9500, and 9501. If you are using the Windows firewall function, set it so that these ports are not blocked.

Important: Network settings

The settings listed here are only an example. The procedure and settings given in this manual do not match every environment. In an environment where the computer being used is connected to multiple network devices and the Internet, inappropriate settings will have a large effect on the entire network. For detailed information about the settings, consult your network administrator.

Determine the IP address, subnet mask, and default gateway to connect the computer and the machine via LAN.

Connecting the Computer and the Machine Directly

	Computer	Machine side	Remarks
IP address	192.168.1.XXX	192.168.1.YYY	Ensure that these are unique numbers compared to other network devices and that "XXX" and "YYY" represent unique numbers from 1 to 254.
Subnet mask	255.255.255.0	255.255.255.0	Set the computer and the machine to the same number.
Default gateway	192.168.1.AAA	192.168.1.BBB	"AAA" and "BBB" represent unique numbers from 1 to 254. This setting may not be necessary. For details, check with your network administrator.

Connecting to an Existing LAN Environment

For each setting value, check with your network administrator.

	Computer	Machine side
IP address	Check with your network administrator.	Check with your network administrator.
Subnet mask	Check with your network administrator.	Check with your network administrator.
Default gateway	Check with your network administrator.	Check with your network administrator.

Step 6 : Setting the Computer's Network

* This step is not necessary when connecting to an existing LAN environment. Proceed to P. 54 "Step 7 : Setting the Machine's Network".

Procedure

- 1 Log on to Windows as the computer's administrator or as a member of the "Administrators" group.
- 2 Display the [Internet Protocol Version 4 (TCP/IP) Properties] window on the computer.

Windows 10

- 1 On the [Start] menu, click [Settings].
- 2 Click [Network & Internet] then [Network and Sharing Center].
- 3 Click [Change adapter settings] from the list on the left side of the window.
- 4 Right-click [Ethernet].
- 5 Click [Properties].

Windows 8.1


- 1 On the [Start] menu, click [PC settings] then [Control Panel].
- 2 Click [Network & Internet] then [Network and Sharing Center].
- 3 Click [Change adapter settings] from the list on the left side of the window.
- 4 Right-click [Ethernet].
- 5 Click [Properties].

Windows 7

- 1 On the [Start] menu, click [Control Panel].
- 2 Click [Network & Internet] then [Network and Sharing Center].
- 3 Click [Change adapter settings] from the list on the left side of the window.
- 4 Right-click [Local Area Connection].
- 5 Click [Properties].

3 Click [Properties].


If the [User Account Control] window appears, click [Continue].

* With a Wi-Fi connection, the name of this window is [Wi-Fi Status].

* On Windows 7, the name of this window is [Local Area Connection Status] or [Wireless Network Connection Status].

4 Select [Internet Protocol Version 4 (TCP/IPv4)], and then click [Properties].


* With a Wi-Fi connection, the name of this window is [Wi-Fi Properties].

* On Windows 7, the name of this window is [Local Area Connection Properties] or [Wireless Network Connection Properties].

2. Connecting to the Computer

- 5 Select [Use the following IP address], and then enter the number you determined in P. 51 “Step 5 : Determining the LAN Address”.

- 6 Click [OK] or [Close] to close the windows that have been displayed in the procedure.

Step 7 : Setting the Machine's Network


Procedure

- 1 VPanel will be displayed.

Click  (the VPanel icon) in the task tray on the desktop.

VPanel will be displayed. If you cannot find  in the task tray, start the program from the Windows [Start] screen (or the [Start] menu).

☞ P. 50 “Step 4 : Selecting the Machine to Connect To”


2 Click .

The [Settings] window is displayed.


3 Click the [Settings] tab, and then click [Add connection].

2. Connecting to the Computer

4 Enter the items.

- Enter a unique name for [LAN-connected printer name]. The name you enter here will be displayed in VPanel and as the name of the computer's printer.
- Enter the values determined in "Step 5 : Determining the LAN Address" (P. 51) for the [IP address], [Subnet mask], and [Default gateway].

Add LAN connection

USB-Connected printer name
DGSHAPE DWX-42W

LAN-Connected printer name
DGSHAPE DWX-42W [IP]

IP address
192 . 168 . 1 . 254

Subnet mask
255 . 255 . 255 . 0

Default gateway
192 . 168 . 1 . 1

OK Cancel

5 Click [OK].


The machine restarts. The settings are complete when the machine finishes restarting.

MEMO

- The name entered for [LAN-connected printer name] is added to [Devices and Printers] in Control Panel.
- The name entered for [LAN-connected printer name] is displayed on the [Settings] tab in the [Settings] window in VPanel.
- If you want to change the machine's network settings, click [Set IP] on the [Settings] tab in the [Settings] window in VPanel.

Step 8 : Disconnecting the USB Cable

- Remove the USB cable from the cable clamp.

Step 9 : Connecting the LAN Cable

- Prepare a LAN cable (Category 5 or higher) separately.
- We recommend 100BASE-TX for the Ethernet.


Important

- Do not bind the LAN cable with the cable clamp.
- Wire the LAN cable so that it does not touch the USB cable or the power cord.

The LAN cable coming into contact with other cables will result in communication errors due to electrical noise.

Connecting the Computer and the Machine Directly

Connecting through a LAN hub


Important: Do not connect any device other than the machine and a computer to the LAN hub.

Do not connect any device other than the designated machine and a computer to the LAN hub being used here. Connecting any other device will prevent the settings from being configured correctly and may cause problems with the network settings of connected devices.

2. Connecting to the Computer


Connecting the machine and the computer directly

MEMO

A direct connection to the machine is possible when using a crossover cable. In addition, if the computer supports AUTO MDI/ MDI-X, a direct connection to the machine is possible when using a straight cable.

Connecting to an Existing LAN Environment

Step 10 : Selecting the Machine to Connect To


Procedure

1 Show VPanel.

Click  (the VPanel icon) in the task tray on the desktop.

VPanel will be displayed. If you cannot find  in the task tray, start the program from the Windows [Start] screen (or the [Start] menu).

☞ P. 50 "Step 4 : Selecting the Machine to Connect To"


2 Click .

The [Settings] window is displayed.


3 Click the [VPanel] tab, and then click [Select].

2. Connecting to the Computer

- 4 From the [Usable machine list], click the image of the machine whose [LAN-connected printer name] was set in step 4 of “Step 7 : Setting the Machine's Network” (P. 54).

The clicked machine is displayed under [Machine to operate] and is displayed dimly in the [Usable machine list]. If you click the same machine again in the [Usable machine list], the machine disappears from the [Machine to operate] field.

Important: When the machine is connected to multiple computers.

If the same machine is selected as the [Machine to operate] on multiple computers, communication may become unstable. Before selecting the machine that was connected this time, check that it is not displayed as the [Machine to operate] in the VPanel window on other computers.

- 5 Click [OK].

3. Performing the Initial Work

Initial Work

Situations When “Initial Work” Is Required

- After the initial installation
- When the installation site is changed

Flow of the “Initial Work”

- STEP 1: Filling the Coolant Tank (P.61)
- STEP 2: Automatic Correction (P.68)
- STEP 3: Spindle Run-In (Long) (P.73)

Step 1 : Filling the Coolant Tank

⚠ CAUTION


Do not overfill or tilt the coolant tank.

The fluid inlet on the coolant tank is open. Overfilling or tilting will cause fluid to spill out.

⚠ CAUTION

When raising or lowering the coolant tank, hold the position indicated in the following figure.

Failing to do so may result in your fingers being pinched, leading to injury.

Items to Prepare Separately


- **Water (soft or purified water)** (P.29)
- **Additive (ZAW-1000D)** (P.31)
- **Chelating reagent (ZCH-100D)** (P.31)

3. Performing the Initial Work


Filling

1. Remove the coolant tank.


1 Click .


The following window is displayed.

2 Open the bottom cover.


- 3 Pull out the coolant tank.

- 4 Remove the coolant tank.

2. Fill the machine with coolant.


- 1 Remove the tray and lid of the coolant tank.

Lift the lid and the tray together to remove them at the same time.


- 2 Install the coolant filters (2).

Orient the coolant filters so that the end with the protrusion is at the front (the end where the coolant tank cap is not present).

3. Performing the Initial Work

3 Check the tightening of the coolant tank cap.

4 Make the coolant.

To make the coolant, mix the water, additive (ZAW-1000D), and chelating reagent (ZCH-100D) at fixed ratios. Use the included measuring cup for measuring.

① Mix water and additive in the ratio of “95:5” (water:additive).

There is no need to stir the solution.

The capacity of the coolant tank is approximately 5 L. To fill the tank, use 4,750 ml of water and 250 ml of additive.

② Mix the water + additive solution and chelating reagent in the ratio of “400:1” (water + additive solution:chelating reagent).

There is no need to stir the solution.

To fill the tank, use 5,000 ml of water + additive solution and 12.5 ml of chelating reagent.

① Water:additive = 95:5

* * When filling the tank: water = 4,750 ml and additive = 250 ml

② Water + additive solution:chelating reagent = 400:1

* * When filling the tank: water + additive solution = 5,000 ml and chelating reagent = 12.5 ml

⚠ CAUTION

Be sure to use the specified additive.

Additives are effective in reducing coolant deterioration and raising the milling efficiency in order to maintain product performance. Furthermore, if additives are not used, the coolant may generate an unpleasant odor.

⚠ CAUTION

Be sure to use the specified chelating reagent.

The chelating reagent works to make calcium and similar substances contained in the water inactive, making it easier to clean off milling powder that adheres to the machine and the cutting tool after milling.

Refer to the appropriate safety data sheet (SDS) for the chemical substances used in the additive and chelating reagent and the safety related to those substances.

To purchase additives, contact your authorized DGSHAPE Corporation dealer or access our website (<http://www.dgshape.com/>).

Important

Add an additional 12.5 ml of chelating reagent to the coolant if:


- It is easy for milling waste to adhere to the inside of the machine.
- The milling time exceeds 8 hours, regardless of whether milling waste adheres to the machine.

☞ P. 61 "Step 1 : Filling the Coolant Tank"


MEMO

The ratios required to create approximately one liter of coolant are written on the label affixed close to the coolant tank.

5 Fill the coolant tank with the coolant.


6 After filling the tank with coolant, put the tray and lid on the coolant tank.

3. Performing the Initial Work


3. Return the coolant tank to its original position.

- 1 Return the coolant tank to its original position.


- 2 Push the coolant tank toward the back of the machine.

Push the coolant tank to the point where you feel a click.

Point:


Align the ▷ label affixed to the right side of the tank with the ◁ label affixed to the machine.

Important

Move the coolant tank slowly. Forcefully moving the coolant tank may cause the coolant to spray out.





- 3 Close the bottom cover.

This completes the procedure for filling the coolant tank.
Continue on to “Step 2 : Automatic Correction” (P. 68)

Step 2 : Automatic Correction

Required Tools

Included items			
			
Rotary axis correction jig	ATC magazine correction jig/ A screw	Detection pin for correction	Cloth for care and maintenance

Important

Handle the detection pin for correction and the dummy pin separately to be absolutely sure that you do not confuse them.


If the detection pin for correction is used as the dummy pin even once, the detection pin for correction cannot be used to provide proper correction. If for some reason the detection pin for correction is used as the dummy pin, a new detection pin for correction will be necessary. Contact your authorized DGSHAPE Corporation dealer or access our website (<http://www.dgshape.com/>).

1. Clean the detection point for use during correction.

- ① Open the front cover.
- ② Remove all milling burs that have been set.

MEMO When the milling bur is difficult to remove (using the milling bur removal jig)

When the milling bur is difficult to remove, you can remove it easily by using the milling bur removal jig.


Insert the tip of the milling bur into the tip (the narrow end) of the milling bur removal jig, and then push the milling bur up.

⚠ CAUTION

Caution: Be careful around the tip and other sharp edges. Doing so may cause injury.


- 3 Use the cloth for care and maintenance to wipe the locations shown in the following figures clean of fluid and dirt.

If any fluid or dirt is present in these locations, correction might not occur properly.

2. Attach the automatic correction jig.

- 1 Attach the ATC magazine correction jig.
- 2 Push the detection pin for correction as far as it will go into stocker number 4 in the ATC magazine.

Important


The tapered (slanted) end is the top of the detection pin for correction. Be careful not to confuse which side is up and which is down.

3. Performing the Initial Work


3 Attach the rotary axis correction jig to hole "2" on the rotary axis.

① Align the recessed portion of the jig with the protrusion on the rotary axis, and then insert the jig.

Ensure that there is no gap between the surfaces. It does not matter which of the two recesses on the jig is aligned to the protrusion.

② Use a hexagonal screwdriver to secure the jig in place with a mounting screw.


4 Close the front cover.

3. Start automatic correction.

MEMO


If the display under [MAINTENANCE] in VPanel is as shown below, ignore it at this point. This indication is handled by P. 73 "Step 3 : Spindle Run-In (Long)" which will be carried out after this procedure.


1 Click .

2 Click the [Maintenance] tab, and then click [Automatic correction].

The following window is displayed automatically.


3 Check that the work displayed on the screen is complete.


3. Performing the Initial Work

4 Click [OK].

Automatic correction starts. The machine status light will flash blue. The remaining work time will be displayed on VPanel.


The automatic correction is finished when the message shown in the figure is displayed.

5 Click [OK].

6 Open the front cover.

7 Remove the correction pin and jig.


8 Close the front cover.

Continue on to "STEP 3: Spindle Run-in (Long)."

Step 3 : Spindle Run-In (Long)

Procedure

- 1 Close the front cover.

MEMO


If a window prompting you to perform collet maintenance is displayed, do so according to the user's manual.

- 2 Click .


MEMO

If the display under [MAINTENANCE] in VPanel is as shown below, click . Clicking this icon displays the same window as clicking [Spindle run-in (long)] in step [3](#). Proceed directly to step [4](#).


3. Performing the Initial Work

- 3 Click the [Maintenance] tab, and then click [Spindle run-in (long)].


The following window is displayed.

- 4 Remove the workpiece, rotary axis correction jig, and ATC magazine correction jig if they are attached.
- 5 Click [OK].
- 6 Open the front cover.


The following window is displayed.


7 Spin the spindle around 10 times in either direction by hand.

8 Close the front cover.

Spindle run-in will begin. The machine status light will flash blue. The remaining work time will be displayed on VPanel.

When the message in the figure is displayed, click [OK].


3. Appendix

Specifications.....	78
Power Rating and Serial Number Location	78
Coordinates	78
Dimensional Drawings.....	79
Specifications	81

Specifications

Power Rating and Serial Number Location


Serial number

The serial number is required for maintenance, servicing, and support. Never peel off the label.

Power rating

Use an electrical outlet that meets the ratings for voltage, frequency, and current given here.


Coordinates

Dimensional Drawings


Machine Dimensions

Unit: mm (in.)


Detection Pin for Correction/Dummy Pin Dimensions

Unit: mm (in.)

Workpiece Mounting Part Dimensions

Unit: mm (in.)

Specifications

		DWX-42W	
Millable material*	Glass ceramics, composite resins		
	Pin type	Width × depth × height: max. 40 × 20 × 20 mm (1.5 × 0.8 × 0.8 in.)	
Operating speed	XYZ: 6 to 1800 mm/min (0.24 to 70.9 in./min.)		
Spindle speed	15,000 to 60,000 rpm		
Rotary axis travel angle	A: ±360°		
Number of milling burs housed	6		
Attachable milling bur	Shank diameter: 3 mm, dedicated milling bur		
Compatible compressed air	0.18 to 0.22 MPa		
Interface	USB		
	Ethernet (10 BASE-T/100 BASE-TX automatic switching)		
Control command sets	RML-1, NC code		
Power requirements	AC 100 to 240 V ± 10%, 50/60 Hz (overvoltage category: II, IEC 60664-1), 2.1A		
Power consumption	Approx. 190 W		
Operating noise	During operation: 70 db (A) or less (when not cutting)		
	During standby: 48 db (A) or less		
External Dimensions	Width × depth × height: 482 × 560 × 600 mm (19.0 × 22.0 × 23.6 in.)		
Weight	61 kg (135 lb.)		
Installation environment	Indoor use at altitudes up to 2000 m Temperature: 5 to 28°C (41 to 82°F) Humidity: 35 to 80%RH (no condensation) Ambient pollution degree: 2 (as specified by IEC 60664-1) Short-term temporary overvoltage: 1440 V Long-term temporary overvoltage: 490 V		
Included items	Power cord, USB cable, manuals, software CD, automatic correction jig, milling burs, hexagonal screwdriver, etc.		

* Even when made from the described materials, some workpieces cannot be milled due to specifications or physical properties. For more information, contact your authorized Roland DG Corporation dealer.

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For EU Countries

This product must be disposed of separately at your local waste recycling center. Do not dispose of in household waste bin.

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